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Guidance and manufacturer's declaration-electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The product is suitable for use in all establishments other than domestic and those directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration-electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge(ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	The device power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration-electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Radiated RF IEC 61000-4-3	3 V/m 80MHz to 2,7 GHz	3 V/m 80MHz to 2,7 GHz	Recommended separation distance: $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P} \quad 80\text{MHz to } 800 \text{ MHz}$ $d = 2,3 \sqrt{P} \quad 800\text{MHz to } 2,5 \text{ GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distance between portable and mobile RF communications equipment and the Device			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz d = $1,2\sqrt{P}$	80 MHz to 800 MHz d = $1,2\sqrt{P}$	800 MHz to 2,7 GHz d = $2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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

敘述	顯示器
按 [PRE-TARE] 鍵。 從最左側的數字開始，輸入預扣重量值	 The image shows a scale's display screen. At the top, there is a 'Back' button and the text 'Pre-tare'. The main display area shows '0.000' followed by 'kg'. Below the display is a 'CLEAR' button. To the right of the display is a numeric keypad with buttons for digits 1-9, 0, and a left arrow, along with an 'Enter' button.
用數字鍵 0-9 輸入預扣重值 例：預扣值 5.0 kg，請按 0-5-0-0-0 例：預扣值 13.5 kg，請按 1-3-5-0-0 按 [Enter] 鍵確認預扣重值	 The image shows two states of the scale's display. The top state shows '0-.---' followed by 'kg', with a 'CLEAR' button below the display. The bottom state shows '5.000' followed by 'kg', with a 'CLEAR' button below the display. Both states have a 'Back' button and 'Pre-tare' text at the top, and a numeric keypad on the right.

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The screenshot shows the main interface of a medical scale. At the top, there are fields for 'Patient:' and 'Nurse:'. Below these, it says 'Awaiting transfer:' with a 'PRE-TARE' button. The central display shows a large '-5.000 kg'. To the right is a 'SCAN' button. Below the weight, there are input fields for 'Height (cm)' and 'BMI'. At the very bottom, it displays 'Max 20 kg', 'Min 0.1 kg', and 'e=0.005 kg'.

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The screenshot shows the 'MEMORY' screen. At the top, there is a 'Back' button. The display shows '0-.-.- kg'. Below the display is a 'CLEAR' button. To the right is a numeric keypad with buttons for 1, 2, 3, 4, 5, 6, 7, 8, 9, 0, and an 'Enter' button. The number '3' is highlighted on the keypad.

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The screenshot shows the 'MEMORY' screen. At the top, there is a 'Back' button. The display shows '0.500 kg'. Below the display is a 'CLEAR' button. To the right is a numeric keypad with buttons for 1, 2, 3, 4, 5, 6, 7, 8, 9, 0, and an 'Enter' button. The 'Enter' button is highlighted.

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The screenshot shows a digital scale interface. At the top, there are fields for 'Patient:' and 'Nurse:'. Below these, it says 'Awaiting transfer:' followed by a 'PRE-TARE' button and a red box containing the number '3'. To the right is a 'SCAN' button. The large central display shows '-0.500kg'. Below this, there are input fields for 'Height (cm)' and 'BMI'. At the bottom, it displays 'Max 20 kg', 'Min 0.1 kg', and 'e=0.005 kg'.

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


This is an identical screenshot to the one above, showing the same digital scale interface with the weight '-0.500kg' and the 'PRE-TARE 3' button highlighted.

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


EXIT General			
Auto Off Time	180s	G-Compensation	-----
Backlight	Mid	URL Host	-----
Buzzer	On	Auto Hold	On
Optional Settings	-----	Auto Transfer	Off
Date/Time	-----		
Wifi Setting	-----		

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Copyright© Charder Electronic Co., Ltd. 版權所有。