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### HM200U

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| <b>Guidance and manufacturer's declaration-electromagnetic emissions</b>                                                                                                                     |                   |                                                                                                                                                                                                   |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>The product is intended for use in the electromagnetic environment specified below.<br/>The customer or the user of the product should assure that it is used in such an environment.</p> |                   |                                                                                                                                                                                                   |
| <b>Emission test</b>                                                                                                                                                                         | <b>Compliance</b> | <b>Electromagnetic environment-guidance</b>                                                                                                                                                       |
| RF emissions<br>CISPR 11                                                                                                                                                                     | Group 1           | The product is suitable for use in all establishments other than domestic and those directly connected to a low voltage power supply network which supplies buildings used for domestic purposes. |
| RF emissions<br>CISPR 11                                                                                                                                                                     | Class A           |                                                                                                                                                                                                   |
| Harmonic emissions<br>IEC 61000-3-2                                                                                                                                                          | Class A           |                                                                                                                                                                                                   |
| Voltage fluctuations<br>/flicker emissions<br>IEC 61000-3-3                                                                                                                                  | Compliance        |                                                                                                                                                                                                   |

## Guidance and manufacturer's declaration-electromagnetic immunity

The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such an environment.

| Immunity test                                                                                       | IEC 60601 test level                                                                                    | Compliance level                                                                                        | Electromagnetic environment-guidance                                                                                                                                                                                                                                           |
|-----------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge(ESD)<br>IEC 61000-4-2                                                       | ±8 kV contact<br>±2 kV, ±4 kV,<br>±8 kV, ±15 kV air                                                     | ±8 kV contact<br>±2 kV, ±4 kV, ±<br>8 kV, ±15 kV air                                                    | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%                                                                                                                                   |
| Electrical fast transient/burst<br>IEC 61000-4-4                                                    | ± 2kV for power supply lines                                                                            | ± 2kV for power supply lines                                                                            | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                            |
| Surge<br>IEC 61000-4-5                                                                              | ±1kV line(s) to line(s)<br>± 2kV line(s) to earth                                                       | ± 1kV line(s) to line(s)<br>± 2kV line(s) to earth                                                      | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                            |
| Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | 0% UT for 0,5 cycle<br>0% UT for 1 cycle<br><br>70% UT(30% dip in UT) for 25cycles<br><br>0% UT for 5 s | 0% UT for 0,5 cycle<br>0% UT for 1 cycle<br><br>70% UT(30% dip in UT) for 25cycles<br><br>0% UT for 5 s | Mains power quality should be that of a typical commercial or hospital environment. If the user of the product requires continued operation during power mains interruptions, it is recommended that the product be powered from an uninterruptible power supply or a battery. |
| Power frequency(50, 60 Hz) magnetic field<br>IEC 61000-4-8                                          | 30 A/m                                                                                                  | 30 A/m                                                                                                  | The product power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                          |

NOTE UT is the a.c. mains voltage prior to application of the test level.

## Guidance and manufacturer's declaration-electromagnetic immunity

The product is intended for use in the electromagnetic environment specified below.

The customer or the user of the product should assure that is used in such and environment.

| Immunity test                 | IEC 60601 test level                                                                                | Compliance level                                                                                    | Electromagnetic environment-guidance                                                                                                                                                                                                                                                                                                |
|-------------------------------|-----------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Conducted RF<br>IEC 61000-4-6 | 3 Vrms<br>150 KHz to 80 MHz<br><br>6 V in ISM bands between 0,15 MHz and 80 MHz<br>80 % AM at 1 kHz | 3 Vrms<br>150 KHz to 80 MHz<br><br>6 V in ISM bands between 0,15 MHz and 80 MHz<br>80 % AM at 1 kHz | Portable and mobile RF communications equipment should be used no closer to any part of the product including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.                                                                                         |
| Radiated RF<br>IEC 61000-4-3  | 3 V/m<br>80MHz to 2,7 GHz                                                                           | 3 V/m<br>80MHz to 2,7 GHz                                                                           | <b>Recommended separation distance:</b><br>$d = 1,2 \sqrt{P}$<br>$d = 1,2 \sqrt{P}$ 80MHz to 800 MHz<br>$d = 2,3 \sqrt{P}$ 800MHz to 2,7GHz<br>Where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in meters (m). |
|                               |                                                                                                     |                                                                                                     | Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup>                                                                                                                                            |
|                               |                                                                                                     |                                                                                                     | Interference may occur in the vicinity of equipment marked with the following symbol:                                                                                                                                                                                                                                               |
|                               |                                                                                                     |                                                                                                     |                                                                                                                                                                                                                                                  |

NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the product is used exceeds the applicable RF compliance level above, the product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the product.

- b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

**Recommended separation distance between  
portable and mobile RF communications equipment and the product**

The product is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the product can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the product as recommended below, according to the maximum output power of the communications equipment.

| Rated<br>maximum<br>output power<br>of transmitter<br><br>W | Separation distance according to frequency<br>of transmitter m |                                               |                                                        |
|-------------------------------------------------------------|----------------------------------------------------------------|-----------------------------------------------|--------------------------------------------------------|
|                                                             | 150 kHz to 80<br>MHz<br><br>$d = 1,2\sqrt{P}$                  | 80 MHz to 800<br>MHz<br><br>$d = 1,2\sqrt{P}$ | <u>800 MHz to 2,7<br/>GHz</u><br><br>$d = 2,3\sqrt{P}$ |
| 0,01                                                        | 0,12                                                           | 0,12                                          | 0,23                                                   |
| 0,1                                                         | 0,38                                                           | 0,38                                          | 0,73                                                   |
| 1                                                           | 1,2                                                            | 1,2                                           | 2,3                                                    |
| 10                                                          | 3,8                                                            | 3,8                                           | 7,3                                                    |
| 100                                                         | 12                                                             | 12                                            | 23                                                     |

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $p$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Copyright© Charder Electronic Co., Ltd.版權所有。