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**MS3500**

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官網：[www.chardermedical.com](http://www.chardermedical.com)

E-mail：[info\\_cec@charder.com.tw](mailto:info_cec@charder.com.tw)

Copyright© Charder Electronic Co., Ltd. All rights reserved.

This user manual is protected by international copyright law. All content is licensed, and usage is subject to written authorization from Charder Electronic Co., Ltd. (hereinafter Charder) Charder is not liable for any damage caused by a failure to adhere to requirements stated in this manual. Charder reserves the right to correct misprints in the manual without prior notice, and modify the exterior of the device for quality purposes without customer consent.

Charder Electronic Co., Ltd.  
No. 103, Guozhong Rd., Dali Dist.,  
Taichung City, 41262 Taiwan

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| Guidance and manufacturer' s declaration-electromagnetic emissions                                                                                                                           |            |                                                                                                                                                                                                                         |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The MS3500 Infant Scale is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment. |            |                                                                                                                                                                                                                         |
| Emission test                                                                                                                                                                                | Compliance | Electromagnetic environment-guidance                                                                                                                                                                                    |
| RF emissions CISPR 11                                                                                                                                                                        | Group 1    | The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.                                         |
| RF emissions CISPR 11                                                                                                                                                                        | Class B    | The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2                                                                                                                                                             | Class A    |                                                                                                                                                                                                                         |
| Voltage fluctuations /flicker emissions IEC 61000-3-3                                                                                                                                        | Compliance |                                                                                                                                                                                                                         |

### Guidance and manufacturer's declaration-electromagnetic immunity

The MS3500 Infant Scale is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

| Immunity test                                                                                       | IEC 60601 test level                                                                                                                 | Compliance level                                                                                                                     | Electromagnetic environment-guidance                                                                                                                                                                                                                                         |
|-----------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge(ESD)<br>IEC 61000-4-2                                                       | <u>± 8 kV contact</u><br><u>± 2 kV, ± 4 kV,</u><br><u>± 8 kV, ± 15 kV</u><br><u>air</u>                                              | <u>± 8 kV contact</u><br><u>± 2 kV, ± 4 kV, ±</u><br><u>8 kV, ± 15 kV air</u>                                                        | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%                                                                                                                                 |
| Electrical fast transient/burst IEC 61000-4-4                                                       | ± 2kV for power supply lines<br>+ 1kV for input/output lines                                                                         | + 2kV for power supply lines<br>+ 1kV for input/output lines                                                                         | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                          |
| Surge IEC 61000-4-5                                                                                 | ± 1kV line(s) to line(s)<br>± 2kV line(s) to earth                                                                                   | + 1kV line(s) to line(s)<br>+ 2kV line(s) to earth                                                                                   | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                          |
| Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | <u>0% UT for 0,5 cycle</u><br><u>0% UT for 1 cycle</u><br><br><u>70% UT(30% dip in UT) for 25 cycles</u><br><br><u>0% UT for 5 s</u> | <u>0% UT for 0,5 cycle</u><br><u>0% UT for 1 cycle</u><br><br><u>70% UT(30% dip in UT) for 25 cycles</u><br><br><u>0% UT for 5 s</u> | Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery. |
| Power frequency(50/60 Hz) magnetic field IEC 61000-4-8                                              | <u>30 A/m</u>                                                                                                                        | <u>30 A/m</u>                                                                                                                        | The device power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                         |

NOTE UT is the a.c. mains voltage prior to application of the test level.

## Guidance and manufacturer' s declaration-electromagnetic immunity

The MS3500 Infant Scale is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that is used in such an environment.

| Immunity test                 | IEC 60601 test level                                                                                            | Compliance level                                                                                                       | Electromagnetic environment-guidance                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Conducted RF<br>IEC 61000-4-6 | 3 Vrms<br>150 KHz to 80 MHz                                                                                     | 3 Vrms<br>150 KHz to 80 MHz                                                                                            | Portable and mobile RF communications equipment should be used no closer to any part of the device including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br><br><b>Recommended separation distance:</b><br>$d = 1,2 \sqrt{P}$<br>$d = 1,2 \sqrt{P}$ 80MHz to 800 MHz<br>$d = 2,3 \sqrt{P}$ 800MHz to 2,5 GHz<br><br>Where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in metres (m).<br><br>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> , should be less than the compliance level in each frequency range <sup>b</sup> .<br><br>Interference may occur in the vicinity of equipment marked with the following symbol:<br><br> |
| Radiated RF IEC 61000-4-3     | <u>6 V in ISM bands between 0,15 MHz and 80 MHz</u><br><u>80 % AM at 1 kHz</u><br><br>3 V/m<br>80MHz to 2,7 GHz | <u>6 V in ISM bands between 0,15 MHz and 80 MHz</u><br><u>80 % AM at 1 kHz</u><br><br>3 V/m<br><u>80MHz to 2,7 GHz</u> |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |

NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.
- b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

### Recommended separation distance between portable and mobile RF communications equipment and the MS3500 Infant Scale

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter<br><br>W | Separation distance according to frequency of transmitter m |                                        |                                         |
|----------------------------------------------------|-------------------------------------------------------------|----------------------------------------|-----------------------------------------|
|                                                    | 150 kHz to 80 MHz<br>$d = 1,2\sqrt{P}$                      | 80 MHz to 800 MHz<br>$d = 1,2\sqrt{P}$ | 800 MHz to 2,5 GHz<br>$d = 2,3\sqrt{P}$ |
| 0,01                                               | 0,12                                                        | 0,12                                   | 0,23                                    |
| 0,1                                                | 0,38                                                        | 0,38                                   | 0,73                                    |
| 1                                                  | 1,2                                                         | 1,2                                    | 2,3                                     |
| 10                                                 | 3,8                                                         | 3,8                                    | 7,3                                     |
| 100                                                | 12                                                          | 12                                     | 23                                      |

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $p$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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