

# MS7800


## Patient Transfer Scale



### Weigh immobile patients quickly and easily

The MS7800 Patient Transfer Scale allows medical personnel to conveniently and reliably measure the weight of patients during the transfer process, saving valuable time and effort. Portable and easy-to-operate, handles can be used for wall storage, saving valuable space, the Patient Transfer Scale is a necessity in every emergency room and ambulance!





<b>Model</b>	<b>MS7800</b>
<b>Capacity</b>	250kg
<b>Graduation</b>	500g
<b>Device Weight</b>	11 kg
<b>Unit</b>	kg or lb
<b>Key Functions</b>	ON/OFF/ZERO/TARE , HOLD/UNITS
<b>Screen Size</b>	28 x 75mm LCD
<b>Power Source</b>	Embedded rechargeable battery
<b>Device Dimensions</b>	1800mm(L) x 700mm(W) x 30mm(H)
<b>Platform Size</b>	1800mm x 510mm



More Info



Measurement of patient weight is essential for administration of correct drug dosages and treatment. Take the guesswork out of the equation, and improve patient safety.

Patients with limited mobility can be difficult to measure. Until now, options have been scales that are undignified for the patient, take too much time to setup and use, or are prohibitively expensive.

We've combined a patient transfer board with a OIML Class III-approved weighing scale, handles can be used for wall storage, saving valuable space. Providing emergency rooms with the perfect space and time-saving solution!

## Features

### Save time

- Measure weight as part of existing patient transfer procedure
- Receive weight results quickly and easily
- Light and portable


### Reduce exertion

- Weigh patients safely, with less physical effort
- Results displayed on large LCD screen
- Non-irritant surface is safe and easily cleaned
- No need to constantly swap batteries - keep the scale ready-to-go using magnetic charging port

### Precise

- 250 kg capacity, 500 g graduation
- Tried and tested in 30 UK hospitals
- Simplify workflow and improve efficiency
- Hold function stabilizes weight of unstable patients


## How is the Patient Transfer Scale used?




**Step 1** Place patient transfer board between trolley and bed (or between two beds)



**Step 2** Assist patient in lying upon transfer sheet and transfer board



**Step 3** Pause to allow scale to measure weight. By using hold function, medical staff can view results afterwards



**Step 4** Complete transfer procedure and remove Patient Transfer Scale

