

MS3830

Wheelchair Scale

The high-capacity MS3830 wheelchair scale is portable and easily transported on rolling castors. Optional handrails can be attached for added stability when used by patients who have difficulty standing. By deducting the weight of the wheelchair using the scale's tare function, a patient's exact weight can be accurately measured.



Image contains optional handrail

Functions

- BMI
- Tare
- Portable Handle
- Print
- Pre-Tare
- Wheel Castors
- Hold
- Backlight
- Data Transfer

Specifications

Indicator	DP3710		
Capacity	400 kg	400 kg	300 kg
Graduation	100 g	200 g	100 g
CE NAWI	—	III	III
Overall Dimensions	1150(W) x 800(D) x 66(H) mm (not including handrails)		
Platform Dimensions	900(W) x 740(D) mm		
Device Weight	32 kg		
Power Supply	AA batteries*6 / Power Adapter / Rechargeable battery pack (optional)		
Optional Accessories	Thermal Printer, Display Stand, Handrails		



Castor wheels for portability



Push wheelchair onto platform via two-way ramps



Quick lock and release for handrails