

CUPID 1 Infant Scale

The CUPID 1 digital infant scale features sloped edges to safely and comfortably hold infants during weight measurement. Highly accurate, the scale can track even slight changes in infants' weight.

The scale is light and portable, well-suited for pediatrician's offices, clinics, hospitals, and primary care centers. It can be used to weigh toddlers as well after removing the measurement tray, and features a hold function to give accurate readings even for moving subjects.





Functions



- Removable Tray



- TARE
- Unit switch (kg/lb)
- HOLD

Model	CUPID I Infant Scale
Capacity	20 kg
Graduation	10g
Overall Dimensions	320(D) x 520(W) x 80(H) mm
Tray Dimensions	320(D) x 520(W) mm
Platform Dimensions	300(D) x 300(W) mm
Device Weight	2.2 kg
Power Supply	9V battery or power adapter



Charder Electronic Co., Ltd.

No.103, Guozhong Rd., Dali Dist., Taichung City 41262, Taiwan

Tel 886-4-2406-3766

Fax 886-4-2406-5612

info_cec@charder.com.tw

www.chardermedical.com