


HM100P

Mechanical Infant Stadiometer



Features

- Lightweight and compact
- Foldable for easy storage
- Quick and simple setup
- Designed for easy transport



Easy-to-read measurement window



Easy Transport for Field Use



AR-00001 Carry Bag (Optional)

Specifications

Measuring range	10-100 cm
Graduation	1mm
Dimensions	1120(W) x 310(D) x 120(H) mm (Open) 63(W) x 310(D) x 120(H) mm (Folded)
Material	ABS Plastic
Weight	1.65 kg

