



Height Measurement

# HMI 10M

## Infant Measuring Mat

Measuring baby length is as easy as laying them down on top of the HMI 10M Infant Measuring Mat. The soft, flexible, vinyl mat is printed on both edges with measurements available in centimeters and inches. The rigid footboard folds to the bottom of the infant's feet, and works with the headboard to ensure accurate measurements, with instructions printed directly on the mat.

Professionals will appreciate that the Infant Measuring Mat is portable, lightweight, and easily stored.



### Specifications

Measuring range 10 - 110 cm (3.9 - 43.3 in)

Graduation 5mm, 1/4 in