

HM200D


Digital Stadiometer

The HM200D digital stadiometer is designed to be wall-mounted, or used with a variety of compatible Charder column scales, featuring a foldable head piece lever with a lock button, and LCD screen with hold function.

High-quality aluminum construction and durable ABS plastic are built to last, providing medical-grade accuracy.

Design Features

- Independent LCD screen with hold function
- Foldable head piece lever with lock button.



Foldable head stopper



Functions

- Hold lock function
- Unit switch (cm/in)
- Folding headpiece

Specifications

Capacity	120-200 cm (47.2-78.7 inch)
Graduation	1 mm (1/16 in)
Overall Dimensions	40 (D) x 70 (W) x 1780 (H) mm
Device Weight	0.85 kg
Power Source	AAA batteries X 2
Mounting	Scale-mounted and Wall-mounted

