




## HM200P Portable Multi-Use Height Stadiometer

The classic portable multi-use Height Stadiometer can be used to measure both vertically and horizontally! Easily broken down into smaller pieces for convenient transportation using the included carry case, it's ideal for traveling nurses, clinics, hospitals, and anywhere that "measurement on the go" is needed!





## FEATURES

- Vertical and horizontal use - just swap out head stoppers!
- Height markers printed on each side for easy viewing
- Convenient and easy to setup and dismantle using carry box
- Large floor plate ensures stability
- Wall stabilizers provide maximum stability
- Light weight and detachable parts make installation and measurement simple



- Can be used to measure infants using optional horizontal mode head stoppers



- Portable carry box



## SPECIFICATIONS

<b>Model</b>	<b>HM200P</b>
<b>Measuring Range</b>	14 – 205 cm / 5.5 – 81 inch
<b>Graduation</b>	1 mm, 1/16 inch
<b>Units</b>	mm, inch
<b>Standard Accessory</b>	Support plate x 1 / Head stopper x 1 / Base platform x 1 / Height rod x 4 / (Optional) Horizontal head stopper x 2 / Horizontal extension rod x 1
<b>Weight</b>	2.0 kg / 4.4 lb (Carry Case 0.6 kg / 1.3 lb)
<b>Dimensions</b>	350(W) x 400(D) x 2430(H) mm Carry Case: 645(W) x 365(D) x 125(H) mm