



# HM210D

## Digital Wall-Mounted Stadiometer

The HM210D professional stadiometer features medical-grade accuracy with a digital LCD screen with manual markings on the height rod. It mounts securely to a wall, and is equipped with a heel plate for accuracy. The headpiece slides easily up and down, and folds flat when not in use.

Switch between metric and imperial units with the click of a button, for international use.



Clear LCD screen



Foldable head stopper



Heel kick plate

### Functions

- Hold lock function
- Unit switch (cm/in)
- Folding headpiece

### Specifications

Capacity	100-210 cm (39-83 in)
Graduation	1 mm (1/16 in)
Overall Dimensions	80(D) x 70(W) x 2100(H) mm
Device Weight	1.25 kg

