



Height Measurement

HM80D

Digital Infant Stadiometer

The HM80D Digital Infant Height Stadiometer easily measures the length of babies as they lie on a scale. Designed to facilitate simultaneous weighing and height measurement, results can be read both using the markings or the embedded LCD screen.



Specifications

Measuring range 35 - 80 cm (14 - 31.5 in)

Graduation 1 mm, 1/16 in

HM80M

Mechanical Infant Stadiometer

Designed to facilitate simultaneous weighing and height measurement, the HM80M can link with compatible Charder infant scales for intuitive and streamlined infant measurement.



Specifications

Measuring range 35 - 80 cm (14 - 31.5 in)

Graduation 1 mm, 1/16 in