

# MG4800

## Grip Strength Dynamometer

The MG4800 medical handgrip dynamometer is specially designed for arm rehabilitation and evaluation. It's a great tool for therapists and rehabilitation trainers to aid in strengthening and training the muscles in the hand and arm.

The dynamometer can be adjusted to different levels of resistance for different individuals. Three modes are available, measuring maximum grip strength, average grip strength, or number of repetitions that meet the selected minimum grip force.



Multiple measurement modes




Multiple spring options

### Functions

- Average grip strength
- Maximum grip strength
- Count Mode

### Specifications

Capacity	80 kg (Blue spring) 40 kg (Green spring) 20 kg (Red spring) 130 kg (Optional)
Grip span	64 mm
Graduation	0.1 kg
Overall Dimensions	212 x 55 x 101 mm
Device Weight	0.33 kg
Power Supply	CR-2450 Battery
Add-ons	Carry case, three spring options



Different capacity of spring



Carrying case