


MHS2510I

Lift Scale



Designed to meet different lift/mobility systems



More joints are available



Features

- Large size LCD display
- Hold function to lock weight of moving patient
- Joint protection design
- Can be used with most manufacturers lifts
- Compact size and reliable performance

Specifications

Model	Display	Capacity	Graduation	Dimensions	Functions	Data Transmission	Joint
MHS2510I	Embedded 1.3 inch LCD screen with 5.5 digits	300kg 400kg	100g 200g	52 x 122 x 180 mm	Power on/off, zero, weight lock, BMI, tare	N/A	SS-6381, SS-6663, SS-7121, SS-7122 more joints are available