




# Versatile Needs

## MHS2510I

### Lift Scale



Designed to meet different lift/mobility systems



More joints are available



## Features

- Large size LCD display
- Hold function to lock weight of moving patient
- Joint protection design
- Can be used with most manufacturers lifts
- Compact size and reliable performance

## Specifications

Model	Display	Capacity	Graduation	Dimensions	Functions	Data Transmission	Joint
MHS2510I	Embedded 1.3 inch LCD screen	300kg 400kg	100g 200g	122(W) x 52(D) x 180(H) mm	Zero/Tare, Hold, BMI	N/A	SS-6381, SS-6663, SS-7121, SS-7122 more joints are available