

MS21NEOV

Digital Column Infant Scale

Features

- Tare function for deduction of weight (ex: diaper)
- Ergonomic tray for comfort and safe weighing
- Tray is non-allergenic and removable for easier cleaning
- Bubble level indicator
- Results can be transferred to PC or thermal printer



Adjustable angle of display



HM80D Digital Infant Height Rod



HM80M Mechanical Infant Height Rod

Functions

- Weight & Height measurement
- BMI
- Hold
- Print
- Tare
- Up to 5 sets of pre-tare memories

Specifications

Indicator	DP3710 Dual LCD
Capacity	15 kg
Graduation	5 g or 2 g < 6 kg > 5 g
Overall Dimensions	560(W) x 450(D) x 470(H) mm
Tray Dimensions	560(W) x 290(D) x 65(H) mm
Device Weight	4.8 kg
Power Supply	AA Batteries*6 / Rechargeable Battery (optional)
Data Transmission	USB or Wireless(optional)

Optional Accessories

Height Measurement	HM80D, HM80M
Thermal Printer	TP2100, TP2110(wireless)
Carrying Bag	AR-00001, AR-2491



AR-00001 carry bag



AR-2491 carry bag