



## MS2 I NEOV Digital Column Infant Scale



Adjustable angle of display



Compatible with Charder infant height stadiometers



AR-2491 Optional carry bag



AR-00001 Optional carry bag

## Features

- Tare function for deduction of weight (ex: diaper)
- Ergonomic tray for comfort and safe weighing
- Tray is non-allergenic and removable for easier cleaning
- Bubble level indicator
- Results can be transferred to PC or thermal printer

### Specifications

### Accessories

Model	Display	Capacity	Graduation	Dimensions	Functions	Data Transmission	Height Rod	Carry Bag	Printer
MS2 I NEOV	DP3710 Dual LCD	15kg	5g or 2g < 6kg > 5g	Overall: 560(W) x 450(D) x 470 (H) mm  Tray: 560(W) x 290(D) x 65(H) mm	Weight & Height measurement, Hold, Tare, BMI, Print  Up to 5 sets of pre-tare memories	USB or Wireless(optional)	HM80M HM80D	AR-2491 AR-00001	TP2100 TP2110 (wireless)