

Nascent Baby & Toddler

MS2400

Low Profile Portable Baby Scale



Image contains optional Infant Stadiometer

Features

- Two aluminum handles for easy transport
- Slim and portable platform design
- Easy to read & clear indicator display
- Memory function for remembering the last weight on hold
- Auto-Hold function to lock weight of moving infant

Functions

- Weight Measurement
- Weight Memory & Recall
- Weight lock
- tare

Specifications

Capacity	20kg
Graduation	10g
Overall Dimensions	640(W) x 355(D) x 60(H) mm
Tray Dimensions	500(W) x 350(D) x 60(H) mm
Device Weight	2.8 kg
Power Supply	AA batteries × 6 / Power Adapter
Accessories	Infant Stadiometer, TP2100 Thermal Printer

Optional Accessories

Height Measurement	HM80D, HM80M
Carrying Bag	AR-248 I



Compatible with Charder infantometer



Clear 1.1" LCD display



Soft sponge handle for easy transport



AR-249 I Optional carry bag

Charder Electronic Co., Ltd.

No.103, Guozhong Rd., Dali Dist., Taichung City 412023, Taiwan
 TEL +886 4 2406 3766 FAX+886 4 2406 5612
 Email: info_cec@charder.com.tw www.chardermedical.com