



Newborn & Toddlers

MS2400

Low Profile Infant Scale



Compatible with Charder infant height stadiometers



Clear 1.1" LCD display



Soft rubber handle for easy transport



AR-2481 Optional carry bag

Features

- Two aluminum handles for easy transport
- Slim and portable platform design
- Easy-to-read indicator
- Memory function for remembering the last weight on hold
- Auto-Hold function to lock weight of moving infant

Specifications

Optional Accessories

Model	Display	Capacity	Graduation	Dimensions	Platform Size	Functions	Height Stadiometer	Carry Bag
MS2400	Embedded 1.1 inch LCD	20kg	10g	640 x 355 x 60 mm (including head pad) Weight: 2.8 kg	500 x 350 x 60 mm	Weight measurement, Recall, Hold, Tare	HM80M HM80D	AR-2491