



Senior & Disabled

MS2504

Handrail Scale



Dual LCD displays weight, height and BMI simultaneously



Ideal for those who require more support during measurement



Castor wheels for mobile use



ST-8236 Optional power adapter pouch



Features

- Ideal for patients who have difficulty stepping onto or balancing on a floor scale
- Integrated castor wheels for easier relocation and storage
- Sturdy handrails provide support for patients
- Non-slip platform for patient safety
- Able to output data to PC or thermal printer

Specifications

Accessories

| Model | Display | Capacity | Graduation | Dimensions | Platform Dimensions | Functions | Data Transmission | Height Rod | Printer | Adapter Pouch |
|--------|--|----------|------------|---------------------|---------------------|--|---------------------------|----------------------------|-----------------------------|---------------|
| MS2504 | DP3710 Dual LCD display with 1 inch screen | 300kg | 100g | 550 x 550 x 1090 mm | 550 x 550 mm | BMI, Print, Weight & Height measurement, Hold, Tare, 5 sets of Pre-tare memories | USB or Wireless(optional) | HM201M HM201D HM200D | TP2100 TP2110 (wireless) | ST-8236 |