

MS2504

Handrail Scale

The MS2504 medical handrail scale's low-profile, anti-slip platform, handrails, and high capacity load cells makes it an ideal weighing scale for patients who need additional support while standing. The platform is reinforced and wide enough to accommodate a chair for patients who cannot stand.




Optional
Ultrasonic Stadiometer and
Printer

Functions

- BMI
- Hold
- Tare
- Pre-Tare
- Bubble Level
- Adjustable Feet
- Castor Wheels
- Printer

Specifications

Capacity	300 kg
Graduation	0.1 kg
Dimensions	550(W) x 550(D) x 1090(H) mm
Device Weight	20.8 kg
Functions	BMI, Hold, Tare, Pre-Tare, Bubble Level, Adjustable Feet, Castor Wheels, Printer
Data Transmission	USB, Wireless (optional)
Optional Accessories	TP2100 Thermal Printer Manual Mechanical Stadiometer: HM201M Manual Digital Stadiometer: HM201D, HM200D Automatic Ultrasonic Stadiometer: HM200U



Optional printer



Castor wheels for mobile use