

# MS3200

## Portable Medical Scale

The MS3200 offers high capacity and durability while still being light and portable. With an oversized measurement platform, the MS3200 is suitable for larger patients, and can be battery-powered for mobile use, or powered by adapter if used indoors for a longer time.



AR-233 I Carrying Bag

### Functions

- Tare
- Hold
- Unit (kg / lb)

### Specifications

Capacity	300 kg / 660 lb
Graduation	100 g / 0.2 lb
Key Functions	TARE, ON/ZERO/OFF, UNIT, HOLD
Overall Dimensions	315(W) x 350(D) x 80(H) mm
Device Weight	3.6 kg (7.9 lb)
Accessories	AR-233 I Carrying Bag
Data Transfer	RS232
Power supply	9V Alkaline battery and Adaptor

## Charder Electronic Co., Ltd.

No.103, Guozhong Rd., Dali Dist., Taichung City 41262, Taiwan

Tel 886-4-2406-3766 Fax 886-4-2406-5612 info\_cec@charder.com.tw www.chardermedical.com