




## MS3500

### Baby Scale with Removable Tray

The MS3500 infant digital scale features sloped edges to safely and comfortably hold infants during weighing. Highly accurate, the scale is able to track even slight changes in infants' weight. With a removable baby tray and 20kg capacity, it can be used to weigh toddlers as well.

Receive height results using the optional stadiometer!

**View height results on LCD screen** (and transfer results)  
when compatible infant stadiometers (HM80DT, HMI00D) are attached



## Functions

- Hold
- Tare
- Removable Tray

## Specifications

Capacity	20 kg
Graduation	5g < 10 kg > 10g
Overall Dimensions	560(D) x 350(W) x 120(H) mm
Platform Dimensions	325(D) x 310(D) x 32(H) mm
Device Weight	2.6 kg
Accessories	Infant Height Stadiometer AR-00001 Carrying Bag
Data Transfer	None
Software	None



Measure toddlers with measurement tray removed

## Option

Height Measurement	HM80D, HM80M ,
	HM80DT, HMI00D



AR-2491  
Optional carry bag for baby scale.