

MS3510

Infant & Toddler Scale C€ 2460

The lightweight and portable MS3510 Wireless Infant & Toddler Scale is OIML certified, enjoying excellent precision and premium-grade features to streamline the weighing process. The curved weighing tray improves measurement safety, and can be removed to measure toddlers as well! Backlit display allows for reading of results even in dim settings, and other practical quality-of-life features like Milk Intake and Hold all make life easier for practitioners.



Capacity	20 kg	20 kg	50 kg
Graduation	5 g	10 g	20 g < 20 kg > 50 g
C€ NAWI	—	III	III
Weight	1.7 kg (without tray) 2.56 kg (20 kg model) 2.69 kg (50 kg model)		
Dimensions	313(W) x 300(D) x 80(H) mm (without tray) 560(W) x 308(D) x 140(H) mm (Overall)		
Power Supply	4 x AA batteries USB type C (Backup)		
Display Type	LCD with white backlight		
Functions	Hold, Tare, Recall, Memory, Milk Intake, Send, Net, Data Transfer, Backlit LCD, OT 3D G-sensor		
Data Transmission	Wireless (Optional)		

Features

- Milk-intake function helps calculate milk intake with a single button press
- Hold function stabilizes weight for easier measurement of fussy infants
- Tare feature deducts weight of diaper, blanket, and other accessories
- Light and portable for easier transportation
- 3D G-sensor automatically detects if platform is level, for accurate measurement



Milk Intake Function



Clear centimeter and inch height markings



AR-2491
Optional carry bag for baby scale.



Measure toddlers with
measurement tray removed

Charder Electronic Co., Ltd.

No.103, Guozhong Rd., Dali Dist., Taichung City 41262, Taiwan
TEL +886 4 2406 3766 FAX +886 4 2406 5612
Email: info_cec@charder.com.tw www.chardermedical.com