


MS3830

Wheelchair & Stretcher Scale

The high-capacity MS3830 Wheelchair & Stretcher Scale is portable and easily transported on rolling castors, and can also be used as a bariatric scale, especially when used with optional handrails for added stability when used by patients who have difficulty standing independently. By deducting the weight of the wheelchair or stretcher using the scale's tare function, a patient's exact weight can be accurately measured.



The basic configuration includes the indicator. The indicator stand, thermal printer, and display handrail are optional.

Functions

- BMI
- Print
- Hold
- Tare
- Pre-Tare
- Backlight
- Portable Handle
- Wheel Castors
- Data Transfer

Specifications

Indicator	DP3710		
Capacity	400 kg	400 kg	300 kg
Graduation	100 g	200 g	100 g
CE NAWI	—	III	III
Overall Dimensions	1150(W) x 800(D) x 66(H) mm (not including handrails)		
Platform Dimensions	900(W) x 740(D) mm		
Device Weight	32 kg		
Power Supply	AA batteries*6 / Power Adapter / Rechargeable battery pack (optional)		
Optional Accessories	Thermal Printer, Display Stand, Handrails		



Carrying handle for easy transport



Castor wheels for portability



Quick lock and release for handrails