

# MS3830

## Wheelchair Scale



The high-capacity MS3830 wheelchair scale is portable and easily transported on rolling castors. Optional handrails can be attached for added stability when used by patients who have difficulty standing.

By deducting the weight of the wheelchair using the scale's tare function, a patient's exact weight can be accurately measured.



Castor wheels for portability



Quick lock and release for handrails (optional)



Push wheelchair onto platform via two-way ramps

### Features

- BMI
- Print
- Hold
- Tare
- Pre-Tare
- Backlight
- Portable Handle
- Wheel Castors
- Data Transfer

Model	MS3830 Wheelchair Scale
Capacity	300 kg
Graduation	100g
Overall Dimensions	1150(D) x 800(W) x 66(H) mm
Platform Dimensions	900(D) x 740(D)
Device Weight	28.6 kg
Power Supply	AA batteries*6 / Power Adapter / Rechargeable battery pack (optional)
Accessories	Thermal Printer, Display Stand, Handrails