

MS3880

Wireless Heavy Duty Wheelchair Scale

The high-capacity MS3880 wheelchair scale is portable and easily transported on rolling casters. Optional handrails can be attached for added stability when used by patients who have difficulty standing.

Results can be transmitted wirelessly to PC, making the MS3880 an excellent option for integration with EMR systems.



Image contains optional handrail



Castor wheels for portability



Push wheelchair onto platform via two-way ramps



Quick lock and release for handrails



Quickly input patient ID via barcode scanner

Functions

- BMI
- Pre-Tare
- Bubble Level
- Hold
- Low-profile Ramp
- Castor wheels
- Tare
- Color Touchscreen
- Handle



Specifications

Capacity	300 kg
Graduation	0.1 kg
Overall Dimensions	1150(W) x 800(D) x 66(H) mm (not including handrails)
Device Weight	28.6 kg
Accessories	Barcode Scanner SM-2714 Indicator Stand Handrails
Data Transfer	USB, Wireless
Software	Smart Data Manager