



MS49 Series

Medical Digital Column Scale



DP3400 for MS4900



DP4400 for MS4970
Feature with Body Surface
Area (BSA)



Optional castors for mobile use



Features

- Bubble level indication for level checking
- Adjustable rubber feet to adjust level
- The equipped wheels ensure the transportation of the scale is easy and labor-saving
- Able to output data to a PC or thermal printer
- Die-casting base structure offers durable performance

Specifications

Options

Model	Display	Capacity	Graduation	Platform Dimensions	Functions	Data Transmission	Height Rod	Printer
MS4900	DP3400 1.6 inch LCD screen	300kg	100g, 200g, 500g	360 x 310 mm	BMI, weight measurement, weight lock, tare	RS232	HM201M HM201D HM200D	TP2100
MS4970	DP4400 1.4 inch LCD screen		100g		BMI, BSA, weight measurement, weight lock, tare	USB	HM201M	TP2100