

MS4970

Medical Column Scale (Castor Wheels Version)

Waist-level column, large LCD screen, and steel platform with anti-slip mat make the MS4970 an economic weighing scale for busy clinics. An optional height stadiometer can be attached to the column to make easy use of the scale's integrated BMI and BSA functions.




optional Height Stadiometer

Functions

- BMI
- BSA
- Hold
- Tare
- Bubble Level
- Adjustable Feet
- Anti-slip Mat
- Print
- Unit (non-OIML models)

Specifications

Capacity	300 kg / 660 lb
Graduation	100g / 0.2 lb
Overall Dimensions	360(W) x 440(D) x 970(H) mm (without Stadiometer)
Device Weight	7.8 kg
Data Transmission:	USB
Optional Accessories	HM201M Mechanical Height Stadiometer TP2100 Thermal Printer



Unit Switch, Zero, Hold, BMI, Tare, Body Surface Area (BSA) Functions



Castor Wheels and Adjustable Feet



Optional Thermal Printer