

MS4971

Medical Column Scale (Castor Wheels Version)

OIML certified for medical-grade accuracy, the MS4971 features a waist-level column, large LCD screen, and steel platform with anti-slip mat. An optional height stadiometer can be attached to make easy use of the integrated BMI and BSA functions.



optional Height Stadiometer



Unit Switch, Tare, Zero, Hold, BMI, Tare, Body Surface Area (BSA) Functions

Functions

- BMI
- BSA
- Hold
- Tare
- Bubble Level
- Adjustable Feet
- Anti-slip Mat
- Print

Specifications

Capacity	300 kg
Graduation	100 g
Overall Dimensions	360(W) x 440(D) x 970(H) mm (without Stadiometer)
Device Weight	7.8 kg
Data Transmission:	USB
Optional Accessories	HM201M Mechanical Height Stadiometer TP2100 Thermal Printer



Optional castors for mobile use



Optional Thermal Printer