



# Professional Column Scales

## MS497I

### Medical Digital Column Scale



Unit Switch, Tare, Zero, Hold, BMI, Tare, Body Surface Area (BSA) Functions



Optional castors for mobile use



MS497I with optional HM201M Manual Height Stadiometer

## Features

- Bubble level indicator
- Adjustable rubber feet to adjust level
- Castor wheels ensure transportation, storage, and relocation is easy and labor-saving
- Able to output data to PC or thermal printer
- Die-casting base structure offers durable performance

### Specifications

### Accessories

Model	Display	Capacity	Graduation	Platform Dimensions	Functions	Data Transmission	Height Rod	Printer
MS497I	DP4600 1.4 inch LCD screen	300kg	100g	360 x 310 mm	BMI, BSA, Weight measurement, Hold, Tare	USB	HM201M HM200D HM201D	TP2100