

MS497I

Medical Digital Column Scale

OIML certified for medical-grade accuracy, the MS497I features a waist-level column, large LCD screen, and steel platform with anti-slip mat. An optional height stadiometer can be attached to make easy use of the integrated BMI and BSA functions.



Functions

- BMI
- BSA
- Hold
- Tare
- Bubble Level
- Adjustable Feet
- Anti-slip Mat
- Print

Specifications

Capacity	300 kg
Graduation	100 g
Overall Dimensions	360(W) x 480(D) x 1100(H) mm (without Stadiometer)
Device Weight	8.1 kg
Data Transmission:	USB
Optional Accessories	HM201M Mechanical Height Stadiometer TP2100 Thermal Printer



Unit Switch, Tare, Zero, Hold, BMI, Tare, Body Surface Area (BSA) Functions



Optional Printer