

# MS497I

## Medical Digital Column Scale

OIML certified for medical-grade accuracy, the MS497I features a waist-level column, large LCD screen, and steel platform with anti-slip mat. An optional height stadiometer can be attached to make easy use of the integrated BMI and BSA functions.



### Functions

- BMI
- Tare
- Anti-slip Mat
- BSA
- Bubble Level
- Print
- Hold
- Adjustable Feet

### Specifications

|                      |                                                                |
|----------------------|----------------------------------------------------------------|
| Capacity             | 300 kg                                                         |
| Graduation           | 100 g                                                          |
| Overall Dimensions   | 360(W) x 480(D) x 970(H) mm<br>(without Stadiometer)           |
| Device Weight        | 7.8 kg                                                         |
| Data Transmission:   | USB                                                            |
| Optional Accessories | HM201M Mechanical Height Stadiometer<br>TP2100 Thermal Printer |



Unit Switch, Tare, Zero, Hold, BMI, Tare, Body Surface Area (BSA) Functions



Optional castors for mobile use