



## MS4980

### Wireless Column Scale with Ultrasonic Height Measurement

The MS4980 is versatile, low-maintenance, and easily operated. From scanning patient ID via barcode scanner, to wireless data transfer, everything is designed to reduce manual data entry, reducing errors and workload.

View results on the easy-to-read color display, and transfer data to your EMR system automatically.

#### Functions

- BMI
- Pre-Tare
- Adjustable Feet
- Hold
- Color Touchscreen
- Tare
- Bubble Level



#### Specifications

|                    |   |
|--------------------|---|
| Capacity           | 300 kg  |
| Graduation         | 0.1 kg  |
| Overall Dimensions | 480(D) x 360(W) x 1100(H) mm<br>(without Stadiometer)   |
| Device Weight      | 8.2 kg  |
| Accessories        | Barcode Scanner   |
| Height Measurement | HM201M: Manual Mechanical Stadiometer<br>HM200D: Manual Digital Stadiometer<br>HM200U: Automatic Ultrasonic Stadiometer |
| Data Transfer      | USB, Wireless   |
| Software           | Smart Data Manager  |

#### Ultrasonic Height Measurement

|            |           |
|------------|-----------|
| Capacity   | 120-200cm |
| Graduation | 0.1 cm    |



Color touchscreen

Barcode Scanner

Image contains optional Ultrasonic Stadiometer and Barcode Scanner