

MS5410

Standard Chair Scale

The Charder MS5410 digital chair scale is an economical weighing option designed to be pushed into patient rooms, so patients that have difficulty standing do not need to walk to weighing stations.

The MS5410 features push handles for ease of transportation, and four swiveling wheels with brakes for safety during measurement.



Functions


- Body Mass Index
- Tare
- Bubble Level
- Swivel Armrest
- Print
- Pre-Tare
- Castor Wheels
- Hold
- Backlight
- Foldable Footrest

Specifications

Capacity	300 kg
Graduation	100 g
Overall Dimensions	1040 (H) x 950 (D) x 640 (W) mm
Seat Dimensions	Height 560 mm Width 440 mm Depth 360 mm Armrest height: 735 mm
Device Weight	22.4 kg
Accessories	TP2100 Thermal Printer (USB) TP2110 Wireless Thermal Printer (Wireless Module needed) ST-8236 Adapter Carry Bag
Data Transfer	USB, Wireless



Swiveling armrests



Foldable footrests



Wheel castor brakes