

MS5440

Heavy Duty Chair Scale

The MS5440 heavy duty digital chair scale is specially designed for larger/bariatric patients. It is intended to be pushed directly into patient rooms, so patients that have difficulty standing do not need to walk to weighing stations.

The dual LCD screens display weight and BMI simultaneously, and can be mounted at the side or rear of the chair as needed.

The MS5440 features push handles for ease of transportation, and four swiveling wheels with brakes for safety during measurement.



Flexible side/rear display positioning for patient and caregiver viewing.



Lifting armrest



Foldable footrests



Swiveling footrest for easier storage



Wheel castor brakes

Functions

- BMI
- Tare
- Bubble Level
- Print
- Pre-Tare
- Wheel Castors
- Hold
- Backlight
- Data Transfer

Specifications

Capacity	300 kg
Graduation	0.1 kg
Overall Dimensions	920(H) x 1190(D) x 690(W) mm
Seat Dimensions	Height 530 mm Width 565 mm Depth 380 mm Armrest height 700 mm
Device Weight	23 kg
Power Supply	Rechargeable battery pack(optional) / AA batteries*6 / PowerAdapter
Accessories	Thermal Printer
Data Transfer	USB, WiFi (optional), BLE (optional)