

MS5460

Compact Chair Scale

The Charder MS5460 compact chair scale is specially designed for patients that have difficulty standing, allowing for the simple and precise weighing of patients.

The MS5460 features collapsible folding frames for easy transportation, and can be assembled easily by one person without tools.



Swiveling armrests



Foldable footrests



Compact size for delivery

Functions

- BMI
- BSA
- Print
- Hold
- Tare
- Pre-Tare
- Unit

Specifications

Capacity	250 kg / 550 lb
Graduation	100 g (0.22 lb)
Overall Dimensions	590(W) x 1020(D) x 980(H) mm
Seat Dimensions	Height 536 mm Width 460 mm
Device Weight	17.6 kg
Accessories	TP2100 Thermal Printer (USB) ST-8236 Adapter Carry Bag
Data Transfer	USB