

MS546I

Compact Chair Scale

The OIML-certified MS546I compact wheelchair scale is specially designed for patients that have difficulty standing. The remote display can be placed at the side or rear of the chair as needed, allowing for the simple and precise weighing of patients.

The MS546I features collapsible folding frames for easy transportation, and can be assembled easily by one person without tools.



Swiveling armrests



Foldable footrests



Compact size for delivery

Functions

- Body Mass Index
- BSA
- Print
- Hold
- Tare

Specifications

Capacity	250 kg / 550 lb
Graduation	100 g (0.22 lb)
Overall Dimensions	590(W) x 1020(D) x 980(H) mm
Seat Dimensions	Height 536 mm Width 460 mm Depth 480 mm
Device Weight	17.6 kg
Accessories	TP2100 Thermal Printer (USB) ST-8236 Adapter Carry Bag
Data Transfer	USB