



# MS546 I

## Compact Chair Scale



Swivelling armrests



Foldable footrests



Compact size for delivery



Tool-free assembly design

## Features

- Four castor wheels with brakes
- Folding frames collapse for easy transportation
- One-person assembly, tool-free clip fastenings
- Able to output data to PC



### Key Functions

- Hold
- Tare / Zero
- Body Surface Area
- Body Mass Index

## Specifications

Model	Display	Capacity	Graduation	Dimensions	Functions	Data Transmission
MS546 I	DP4600 1.4 inch LCD screen	250kg	100g	960(H) x 1020(D) x 590(W) mm	BMI, BSA, Hold, Weight measurement, Tare	USB