

MS575 I

Portable Medical Floor Scale

The OIML-certified MS575I features medical-grade accuracy, a wide LCD screen, die-cast platform, and anti-slip surface for user safety. Integrated Body Mass Index (BMI) and Body Surface Area (BSA) functions increase functionality for medical clinics.

For traveling clinicians requiring a robust and high-capacity scale, the MS575I is an excellent choice!

Results can be transferred to PC via USB cable for data management.



Functions

- BMI
- Body Surface Area
- Print
- Hold
- Tare
- Bubble Level
- Adjustable Feet
- Anti-slip mat

Specifications

Capacity	300 kg
Graduation	0.1 kg
Overall Dimensions	360(D) x 430(W) x 110(H) mm
Device Weight	5.5 kg
Data Transfer	USB
Accessories (Optional)	AR-278I Carry Bag TP2100 Thermal Printer



Robust handle for easy portability



AR-278I Carry Bag (Optional)