

MS5980

Wireless Infant Height & Weight Scale

Ergonomic and streamlined, the MS5980 is designed to make weighing fast and comfortable for infants. The extra-large tray makes measurement safe even for larger babies, and the integrated height meter automatically sends results to color touchscreen indicator.

Featuring key functions like hold and tare to make it easier to measure babies, a single button press sends results to an EMR system once measurement is complete.



Image contains optional Infant Stadiometer and Barcode Scanner

Functions

- BMI
- Pre-Tare
- Bubble Level
- Hold
- Color Touchscreen
- Adjustable Feet
- Tare
- Wireless Data Transfer

Specifications

Capacity	15kg	20kg (Optional)
Graduation	2g < 6 kg > 5g	5g < 10kg > 10g
Overall Dimensions	670(D) x 470(W) x 360(H) mm	
Device Weight	4.1 kg	
Accessories	Infant Height Stadiometer Barcode Scanner	
Data Transfer	USB, Wireless	



Quickly input subject ID via barcode scanner



Anti-tilt design

Optional Accessories

Height Measurement	HM80D, HM80M
	HM80DT, HM100D