

MS6111

Double-sided Medical Floor Scale

The MS6111 is an affordable, accurate, professional medical scale, designed with a low profile platform to make it easier to step on and off. The double-sided LCD screen makes it possible for results to be viewed easily by either the patient or the operator. When used with the HM202P Portable Stadiometer, weight and height measurement can be conducted simultaneously, saving valuable time and effort.

The MS6111 is also designed to help measure the weight of infants or toddlers who cannot stand on the scale independently! By using the built-in Mother-child function, the mother's weight is automatically subtracted from the total, displaying child weight clearly.



AR-2781 Carry Bag



Functions

- BMI
- Tare
- Print
- Anti-slip Mat
- Hold
- Unit

Specifications

Capacity	250 kg
Graduation	10g < 20kg > 100g
Overall Dimensions	345(W) x 350(D) x 45(H) mm
Device Weight	2.5 kg
Accessories (Optional)	<ul style="list-style-type: none"> • AR-2781 Carry Bag • HM202P Height Stadiometer • TP2110 Thermal Printer



Double-sided indicator for easy viewing



Compatible with HM202P Stadiometer – no cumbersome base attachment needed