

CUPID 3 baby to toddler scale

The multi-function digital baby scale for family



CUPID 3 baby to toddler scale

FEATURES

- Easy-to-read large 1-inch LCD display with blue backlight
- Tare feature that accounts for extra weight attributed to diaper, clothing or blanket
- It includes power saving feature ease your pain of replacing batteries
- Sleek contour design makes baby feel safe and comfortable
- Auto-switch dual direction LCD display once weight up to 20kg



- Milk-intake function helps mother to measure how much milk baby intake with just one step



- Embedded height stopper measure height up to 160cm, 5ft 3inches



- Baby mommy mode is designed for mommy get baby's weight while cuddling the baby.



SPECIFICATIONS

Capacity	0~20 kg / 44 lbs , 20~100 kg / 220 lbs
Graduation	0.01 kg / 0.02 lb , 0.05 kg / 0.1 lb
Measuring unit	kg or lb
Weight	1.6 kg
Dimensions	platform 290x326x41 mm
	tray 600x354x91.5 mm
Power supply	3 x AAA batteries
Display type	LCD with blue backlit