

MHS2710


Lift Scale



Discover the technology of wireless transmission and lithium battery cartridge

Features

- Higher battery capacity up to 2600mAh
- Continuous work for 120 hours
- Lithium battery cartridge design
- Built-in wireless transmission
- Joint protection design



Rechargeable lithium battery cartridge

Specifications

Options

Model	Display	Capacity	Graduation	Dimensions	Functions	Data Transmission	Joint	Printer
MHS2710	Embedded 1.3 inch LCD screen with 5.5 digits	300kg 400kg	100g 200g	60 x 122 x 180 mm	Power on/off, zero, weight lock, BMI, tare	Wireless	SS-6381,SS-6663, SS-7121,SS-7122 more joints are available	TP2110 (wireless)