



MS2 | NEOV Digital Column Baby Scale



Adjustable angle of display



Compatible with Charder infantometer



AR-2491 Optional carry bag



AR-00001 Optional carry bag

Features

- Tare function to tare the weight of diaper and get net weight of baby
- Ergonomic tray shape for baby to snuggle down safely during weighing
- Removable tray is hygienic and easy to clean
- Able to output data to a PC or thermal printer
- Bubble level indication for level checking

Specifications

Options

Model	Display	Capacity	Graduation	Dimensions	Tray Size	Functions	Data Transmission	Height Rod	Carry Bag	Printer
MS2 NEOV	DP3710 Dual LCD display with 1 inch screen	15kg	5g or 2g < 6kg > 5g	560 x 450 x 470 mm	560 x 290 x 65 mm	Weight & height measurement, weight lock, tare, BMI, print out button, M1-5 sets of pre-tare memories	USB or Wireless(optional)	HM80M HM80D	AR-2491 AR-00001	TP2110 (wireless)