

MS345 I Digital Column Scale



TARE
BSA

Body Surface Area(BSA)

The total surface area of the human body.
For many clinical purposes BSA is a better indicator of metabolic mass than body weight because it is less affected by abnormal adipose mass.
BSA is commonly used in medicine, particularly to calculate doses of chemotherapeutic agents and index cardiac output.



Scale with stadiometer HM201M attached



Durable metal platform

Features

- Integrated BSA and BMI function
- Sturdy metal platform
- Low battery consumption
- Non-skid platform surface
- Measures in kilograms or pounds
- Column mounted display and controls
- Able to output data to a PC

Specifications

Options

Model	Display	Capacity	Graduation	Platform Dimensions	Functions	Data Transmission	Height Rod
MS345 I	DP4600 1.4 inch LCD display	300kg	100kg	310 x 310 mm	BMI, BSA, weight measurement, weight lock, tare	USB	HM201M