



## MS38 Series Portable Wheelchair Scale



Handrail Type C



Castor wheels for mobile use



Handrail Type A



Handrail Type B



Quick lock and release system



Foldable to transport



Handrail Type C-Foldable type

## Features

- Low profile platform
- Sturdy handrail supports for geriatric and bariatric patients
- Remote indicator with preset-tare, and semi-auto tare
- Easily accessible for wheelchairs
- Bubble level indication for level checking
- Adjustable rubber feet to adjust level
- Able to output data to a PC or thermal printer

### Specifications

### Options

| Model  | Display                                    | Capacity | Graduation | Dimensions         | Platform Dimensions | Functions                                                                                  | Data Transmission         | Printer           | Adaptor Pouch | Handrail   | Display Stand |
|--------|--------------------------------------------|----------|------------|--------------------|---------------------|--------------------------------------------------------------------------------------------|---------------------------|-------------------|---------------|------------|---------------|
| MS3830 | DP3710 Dual LCD display with 1 inch screen | 300kg    | 100g       | 1150 x 800 x 65 mm | 900 x 740 mm        | BMI, Print, weight & height measurement, weight lock, tare, M1-5 sets of pre-tare memories | USB or Wireless(optional) | TP2110 (wireless) | ST-8236       | SM-3461(R) | SM-2711       |
| MS3850 | DP4400 1.4 inch LCD screen                 |          |            |                    |                     | BMI, weight measurement, weight lock, tare                                                 | USB                       | N/A               |               | SM-3462(L) |               |

Type A: SM-3461(R), SM-3462(L), SS-8444(M) Type B: SM-3461(R) Type C: SM-00001(R)