



MS49 Series

Medical Digital Column Scale



DP3710 for MS4910



DP3400 for MS4900



DP4400 for MS4970
Feature with Body Surface Area (BSA)



Optional castors for mobile use



Features

- Bubble level indication for level checking
- Adjustable rubber feet to adjust level
- The equipped wheels ensure the transportation of the scale is easy and labor-saving
- Able to output data to a PC or thermal printer
- Die-casting base structure offers durable performance

Specifications

Options

Model	Display	Capacity	Graduation	Platform Dimensions	Functions	Data Transmission	Height Rod	Printer
MS4910	DP3710 Dual LCD display with 1 inch screen	300kg	50g, 100g	360 x 310 mm	BMI, Print, weight & height measurement, weight lock, tare, M1-5 sets of pre-tare memories	USB or Wireless(optional)	HM201M HM201D HM200D	TP2110 (wireless)
MS4900	DP3400 1.6 inch LCD screen		100g, 200g, 500g		BMI, weight measurement, weight lock, tare	RS232	HM201M HM201D HM200D	TP2100
MS4970	DP4400 1.4 inch LCD screen		100g		BMI, BSA, weight measurement, weight lock, tare	USB	HM201M	N/A