



MS58 Series Standard Chair Scale



Clear 1.6 inch LCD display



Free swivelling armrests



Foldable footrests

Features

- Revolvable arm-rest allows users easily and comfortably access
- Lockable rear wheels for safety and security
- Able to output data to a PC or thermal printer
- Non slippery seat
- Inclusive of adaptor pouch for mobile use

Specifications

Options

Model	Display	Capacity	Graduation	Dimensions	Functions	Data Transmission	Printer
MS5810	DP3400 1.6 inch LCD screen	200kg	100g	920 x 950 x 630 mm	BMI, print, hold, tare	RS232	TP2100
MS5811		150kg	100g < 100kg > 200g				