



Nascent Baby & Toddler

MS5900

Large Platform Baby Scale



Tray mimics the shape of cradling arms for infant comfort



Anti tilting tray design for optimum safety during weighing



Clear LCD display



ST-358 | Optional carry bag



AR-0000 | Optional carry bag

Features

- Anti-tilt design to prevent infant fall from tray
- Able to output data to a PC or thermal printer
- Weighing memory recall function
- Bubble level indication for level checking
- Adjustable rubber feet to adjust level

Specifications

Options

Model	Display	Capacity	Graduation	Dimensions	Tray Size	Functions	Data Transmission	Height Rod	Carry Bag	Printer
M5900	Embedded 0.9 inch LCD screen with 5 digits	15kg, 20kg	2g < 6kg > 5g, 5g, 5g < 10kg > 10g	670 x 330 x 125 mm	630 x 250 x 70 mm	Weight measurement, weight lock, weight memory & recall, tare Milk In Take	USB or Wireless(optional)	HM80M HM80D	ST-358 AR-0000	TP2110 (wireless)